

**High Commission of Pakistan
Male' (Maldives)

TENDER NOTICE

Bids are invited from the interested parties for supply of the following items at Pakistan High Commission, Male' (Maldives). The minimum specifications and quantity of items required are as follows:

A- Desktop Computer

Description	Quantity
Dell Optiplex 3000 MT Desktop / Intel Core i5-12500 (6 Cores/18MB/12T3.0GHz to 4.6GHz/65w), 16GB DDR4 RAM, 256GB PCIe NVMe SSD, Intel UHD Graphics 630, 1 x Display Ports, 1 x HDMI, Dell Wired USB Keyboard & Mouse	51
Dell S2421HN 24" LED Monitor / Screen size 23.8", Panel Type: IPS, Resolution: FHD (1080p) at 75 Hz.	51
Genuine Windows 11 Pro	51

B- Laptops

Description	Quantity
HP 15-fd0054TU i7-1355U 8GB D4 512GB Intel Iris Xe 15.6" FHD W11	07
Logitech M185 Compact Wireless Mouse	07
HP Prelude Backpack 15.6"	07
Genuine Windows 11 Pro Key	07

C- Television

Description	Quantity
Smart TV 55"	06

D- Computer Table & Chair

Description	Quantity
Computer Table	15
Chair	15

E- Medical Instrument

Sl: #	Description	Quantity
1	K WIRE IMPLANT & INSTRUMENT SET (1 Set)	
	Wire Tightener Handle and Two Pegs	1
	Plier cum Bender cum Cutter	1
	Wire Bending Pliers, (Maun) 6"	1
	Wire Passer, Medium	2
	Vice Grip	1
	Forceps for Holding Cerclage Wire	2
	IMPLANTS	
	Cerclage Wire with Eye of dia. 1.0, length 280 mm	10
	Cerclage Wire with Eye of dia. 1.25, length 280 mm	20
	Suture Wire, each of 1.0 and 1.25 mm dia.	2
	K-Wire Dia. 1.0, Length 150 mm	10
	K-Wire Dia. 1.25, Length 150 mm	10
	K-Wire Dia. 1.5, Length 150 mm	20
	K-Wire Dia. 2.0, Length 150 mm	10
	K-Wire Dia. 2.5, Length 150 mm	10
	K-Wire Dia. 3.0, Length 150 mm	10
	Aluminium Graphic case	1
2	POWER DRILL FOR ORTHOPEDIC SURGERIES	2
	Power source	Battery-powered
	Procedure	Orthopedic surgery
	Sound intensity	70 dB
	Speed	1,200 rpm (7,540 rad.min-1)
	Chuck size	Min.: 0.6 mm (0 in)
		Max.: 8 mm (0.3 in)
	Technical Parameters	
		With reversing function
		Center clamping range:0.6-6mm
		Rotating speed : >=1200rpm
		Input voltage: 110-220V 50/60HZ
		The output voltage: 14.4V
3	ENDERS' NAIL INSTRUMENT SET	
	2.0/2.5/3.0/3.5 MM ENDERS NAIL IMPACTOR DRIVER	1
	4.0/4.5 MM ENDERS NAIL IMPACTOR DRIVER	1
	2.0/2.5/3.0/3.5 MM ENDERS NAIL BENDER	1
	4.0/4.5 MM ENDERS NAIL BENDER	1
	ENDERS NAIL EXTRACTOR HOOK (PAIR)	1
	HOLLOW HAMMER	1

		HANDLE FOR EXTRACTOR HOOK	1
		BONE AWL CURVED	1
		ENDERS NAIL IMPACTOR	1
		INSTRUMENT BOX WITH IMPLANT TRAY	1
4	4.5mm Cannulated Screw Instrument Set		
		T-handle With Quick Coupling	1
		Universal Drill Guide Ø4.5/Ø3.2mm	1
		Screwdriver Shaft With Quick Coupling, SW3.5	1
		Depth Gauge	1
		Cleaning Wire Ø1.6mm	1
		Kirschner Wire Ø1.6mm	10
		Countersink Drill Bit	1
		Parallel Locator Ø1.6mm	1
		Hexagonal Screwdriver SW3.5	1
		6.5 Holding Sleeve	1
		Locator Ø9.5/Ø7.0mm	1
		Drill Bushing Ø7.0/Ø3.2mm	1
		Drill bushing Ø3.2/Ø1.6mm	1
		Trocar Ø1.6mm	1
		Cannulated Drill Bit, Ø3.2mm	1
		Cannulated Tap For Cancellous Bone screwφ4.5mm	1
		Screw Tray	1
		Instrument Tray(CanNUS-4.5)	1
5	Reconstruction Plates, 3.5 mm		20
	these plates are used with 3.5 mm cortex screws;		
	Main application is for use with pelvic and acetabular frx;		
	This may be also used in fractures of the distal humerus (dorsolateral aspect), the clavicle, and the calcaneus;		
	Reconstruction plates have notches alongside the plate, which enables bending in three dimensions;		
	Bending more than 15 deg at any one site should be avoided;		
	Oval holes permit some self-compression if the screw holes are placed eccentrically;		
	Holes accept 4.0 mm cancellous bone screws if the plate is placed over cancellous bone;		
	Screws can be inserted at an angle of approx 25 deg longitudinally and 7 deg sideways;		
	Specification/Measurements		
	Thickness 2.8 mm		
	Width 10 mm		
	Hole spacing 12 mm		
	Note that already low stiffness of the plate is further diminished by bending, & if strong curvature is needed, consider a prebent plate;		
6	Plate Holding Forceps Set		
	Sizes and Lock types (Ratchet)		1
	1.5/2.0 mm		1
	2.4/2.7 mm		1
	3.5 mm		1
	Sizes and Lock types (Spinlock)		1
	1.5/2.0 mm		1
	2.4/2.7 mm		1
	3.5 mm		2
7	Large Fragment Instruments Set		
	Drill Bit (Quick Coupling), 3.2mm x 130mm		1
	Drill Bit (Quick Coupling), 4.5mm x 130mm		1
	Counter Sink for 4.5 & 6.5mm Screws (8.0mm Head)		1
	Quick Coupling Handle (T-Type)		1
	Bone Tap with Quick Coupling for 4.5mm Cortex Screws, L:130mm		1
	Bone Tap with Quick Coupling for 6.5mm Cancellous Screws		1
	Insert Drill Sleeve 4.5/3.2mm		1
	Drill & Tap Sleeve combined 3.2/4.5mm		1
	Drill & Tap Sleeve combined 3.2/6.5mm		1
	Screw Driver Shaft, Q.C. End - 3.5mm Tip		1
	Hexagonal Screw Driver with Sleeve, 3.5mm Tip (Excel)		1
	Depth Gauge for 4.5mm & 6.5mm Screws		1
	Sharp Hook		1
	Large Neutral & Loaded Drill Guide		1
	Pointed Drill Guide		1
8	Small Fragment Locking Plate Instrument Set (AO)		
	Depth Gauge		1
	Guide Pin Sleeve Ø1.5		2
	Drill Sleeve Ø3.2		3
	Small Head Drill Sleeve Ø3.2		3
	Torque Handle 1.5N.m		1
	Quick coupling Hex Screw Driver SW2.5		2
	Lifting Drill		1
	Drill guide Wrench		1
	Self-Holding Screw Driver SW2.5		1
	Double-end Guide Ø3/4		1
	Quick Coupling Drill Bit Ø3.2		3
	Quick Coupling Drill Bit Ø3		2

	Quick Coupling Self-tapping Tap HA4.0	1
	Quick Coupling Locking Screw Tap HC4.0	1
	Quick Coupling Self -holding Screwdriver SW2.5	1
	Quick Coupling Retractor for Slipping Screws	1
	T Quick Coupling Handle	1
	Straight Quick Coupling Handle	1
	Connection for Quick Coupling Handle	1
	Guide Pin Ø1.5	3
	Thread Guide Pin Ø1.5	3
	Postion Stopper Wrench	1
	Screw Holder	1
	Periosteal Dissector 8mm/9mm	1
	Periosteal Dissector 12mm	1
	Dissector	2
	Quick Coupling Countersink Drill	1
	Quick Coupling Hollow Saw For Broken Screw	1
	Sharp End Reduction Forcep	1
	Plate Bender	2
	Self-centering Bone Hold Forcep	2
	Reduction Forcep	1
	Instrument box	1
	Other Surgical Instruments	
9	Description	Quantity
	Balfour Self Retaining Retractor (For use in laparotomy procedures; self-retaining; adjustable instrument).	1
	Richardson Eastman Retractors Double Ended Set of Two 10" (surgical instrument)	1
	Mathieu Rectal Speculum	1
	Percutaneous suprapubic trocar with an outer sheath and inner obturator	2 set
	Circumplast Circumcision Device (Disposable)	50 set

i. **Bid-Security Declaration Form given below is required to be submitted along with respective bid.**

<u>Bid-Security Declaration</u>
<p>The Head of Chancery High Commission of Pakistan, <u>Maldives</u></p> <p>Reference: Tender Notice _____ Our Bid No _____</p> <p>Dear Sir,</p> <p>I / We,..... declare as under:</p> <p>The Bid-Security Declaration is submitted in terms of the above Tender Notice. I/we hereby accept that I/we may be disqualified from bidding for any contract with you for a period of three years from the date of disqualification as may be notified by you, if,</p> <ul style="list-style-type: none"> i. I am/we are in a breach of any of the obligations under the bid conditions. ii. I/we have withdrawn or unilaterally modified/amended/revised, my/our Bid during the bid validity period specified in the form of Bid or extended period, if any. iii. On acceptance of our bid by the High Commission of Pakistan, I/we failed to execute the agreement or fails to provide the goods in accordance with the terms and conditions. <p>Signature:</p> <p>Name & designation Company / Business Name Dated</p>

- ii. Interested parties should submit their bids to the High Commission of Pakistan, R11-02, Kinbigas Magu, Golahi No.1, Hulhumale' Phase-I, Maldives during office hours (0900 – 1600) latest by **19 May 2024 till 1100 hours.**
- iii. All the tenders will be opened in the presence of bidders or their representatives, who may choose to be present at the High Commission of Pakistan, R11-02, Kinbigas Magu, Golahi No.1, Hulhumale' Phase-I, Maldives on **19 May 2024, at 1130 hours.**
- iv. The High Commission reserves the right to change/modify specifications / quantity of the items.
- v. The High Commission also reserves the right to accept or reject any tender as per PPRA Rules.
- vi. The items are required to be delivered at Pakistan High Commission, Hulhumale', Maldives.

(Zulqurnain Ahmad)
Deputy Head of Mission/HOC
Telephone No: 00960 3323005
Email: pahicmale@mofa.gov.pk